

Below is a poorly done tape job. The connector nut is exposed and one of the strain-relief screws has punched thru the tape. Only one layer of tape is used here.



Figure 3 Poorly Done - Single Layer

For more information please contact REF TEK Customer Support at support@reftek.com.

